

Microflex®  
**midknight™**  
 Powder-Free Black  
**Nitrile Exam Gloves**

**Item #**

- Extra Small..... MK-296-XS
- Small ..... MK-296-S
- Medium..... MK-296-M
- Large..... MK-296-L
- Extra Large..... MK-296-XL
- XX-Large..... MK-296-XXL



**Recommended for use in the following industries:**

Automotive / Dental / EMS / Funeral / Healthcare / Industrial Safety / Laboratory / Law & Corrections / Tattoo

**When enhanced tactile sensitivity and consistent grip are a must, make MidKnight™ your glove of choice.**

**MidKnight™** gloves offer the superior combination of durability and protection against day-to-day hazards, in a comfortable nitrile glove. The glove of choice for individuals who demand uncompromising quality and protection, **MidKnight™** gloves provide an exceptional balance of performance, comfort and durability. Designed specifically with a distinct color that masks stains.

**Recommended for use when:**

- Day-to-day hazards require the best combination of durability and consistency
- Working in environments that demand firm grip and reliability
- Precise contact requires enhanced tactile sensitivity

Microflex® MidKnight™ Profile			For current specification information, please visit <a href="http://www.microflex.com">www.microflex.com</a> .	
	ASTM D 6319	MidKnight™		
<b>Length (mm/in):</b>	230 / 9.1	245 / 9.6	<b>Material:</b> Nitrile	<b>Type:</b> Non-sterile
<b>Cuff Thickness (mils):</b>	n/a	3.1	<b>Size:</b> XS - XXL	<b>Color:</b> Black
<b>Palm Thickness (mils):</b>	min. 2.0	4.7	<b>Shape:</b> Ambidextrous	<b>Usage:</b> For single use only
<b>Finger Thickness (mils):</b>	min. 2.0	5.5	<b>Cuff:</b> Standard examination length (9.5"), Beaded	
<b>Tensile Strength:</b>			<b>Glove Exterior:</b> Fully Textured	
Before Aging (MPa):	min. 14	27	<b>Glove Interior:</b> Powder-Free, Chlorinated	
After Aging (MPa):	min. 14	24	<b>Packaging:</b> 1,000 gloves/case, by weight; 10 boxes/case	
<b>Elasticity/Elongation:</b>			100 gloves per dispenser box, by weight	
Before Aging (%):	min. 500	500	<b>Caution:</b> Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.	
After Aging (%):	min. 400	400		

**The reason is the quality... the difference is Microflex®.**

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit [www.microflex.com](http://www.microflex.com). For the name of a distributor in your area, call: **(800) 876-6866**.

